

兔抗 HRH2 多克隆抗体

中文名称：兔抗 HRH2 多克隆抗体

英文名称：Anti-HRH2 rabbit polyclonal antibody

别名：H2R

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HRH2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Histamine is a ubiquitous messenger molecule released from mast cells, enterochromaffin-like cells, and neurons. Its various actions are mediated by histamine receptors H1, H2, H3 and H4. Histamine receptor H2 belongs to the family 1 of G protein-coupled receptors. It is an integral membrane protein and stimulates gastric acid secretion. It also regulates gastrointestinal motility and intestinal secretion and is thought to be involved in regulating cell growth and differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB

Name of antibody:	HRH2
Immunogen:	Synthetic peptide of human HRH2
Full name:	histamine receptor H2
Synonyms :	H2R
SwissProt:	P25021
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	40 kDa
WB Positive control:	Human fetal muscle tissue
WB Recommended dilution:	500-2000

