

兔抗 HAS3 多克隆抗体

中文名称：兔抗 HAS3 多克隆抗体

英文名称： Anti-HAS3 rabbit polyclonal antibody

别名： hyaluronan synthase 3;

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： HAS3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is involved in the synthesis of the unbranched glycosaminoglycan hyaluronan, or hyaluronic acid, which is a major constituent of the extracellular matrix. This gene is a member of the NODC/HAS gene family. Compared to the proteins encoded by other members of this gene family, this protein appears to be more of a regulator of hyaluronan synthesis. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC

Name of antibody:	HAS3
Immunogen:	Fusion protein of human HAS3
Full name:	hyaluronan synthase 3
SwissProt:	O00219
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	63 kDa
WB Positive control:	Human cervical cancer tissue and A549 cell lysates
WB Recommended dilution:	500-2000



