

兔抗 HRK 多克隆抗体

中文名称：兔抗 HRK 多克隆抗体

英文名称：Anti-HRK rabbit polyclonal antibody

别名：DP5; HARAKIRI

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HRK

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the BCL-2 protein family. Members of this family are involved in activating or inhibiting apoptosis. The encoded protein localizes to intracellular membranes. This protein promotes apoptosis by interacting with the apoptotic inhibitors BCL-2 and BCL-X(L) via its BH3 domain. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	HRK
Immunogen:	Synthetic peptide of human HRK
Full name:	harakiri, BCL2 interacting protein (contains only BH3 domain)
Synonyms:	DP5; HARAKIRI

SwissProt:	O00198
IHC positive control :	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	50-200

