

## 兔抗 MSI1 多克隆抗体

中文名称：兔抗 MSI1 多克隆抗体

英文名称：Anti-MSI1 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：MSI1

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein containing two conserved tandem RNA recognition motifs. Similar proteins in other species function as RNA-binding proteins and play central roles in posttranscriptional gene regulation. Expression of this gene has been correlated with the grade of the malignancy and proliferative activity in gliomas and melanomas. A pseudogene for this gene is located on chromosome 11q13.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	MSI1
<b>Immunogen:</b>	Synthetic peptide of humanMSI1
<b>Full name:</b>	musashi RNA-binding protein 1
<b>SwissProt:</b>	O43347
<b>ELISA Recommended diluti</b>	1000-2000

<b>on:</b>	
<b>IHC positive control:</b>	Human liver cancer and Human breast cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	Mouse brain tissue
<b>WB Recommended dilution:</b>	200-1000



