

兔抗 PTPN3 多克隆抗体

中文名称：兔抗 PTPN3 多克隆抗体

英文名称： Anti-PTPN3 rabbit polyclonal antibody

别名： PTPH1; PTP-H1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PTPN3

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This protein contains a C-terminal PTP domain and an N-terminal domain homologous to the band 4.1 superfamily of cytoskeletal-associated proteins.
Applications:	ELISA, IHC
Name of antibody:	PTPN3

Immunogen:	Fusion protein of human PTPN3
Full name:	protein tyrosine phosphatase, non-receptor type 3
Synonyms :	PTPH1; PTP-H1
SwissProt:	P26045
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain and human ovarian cancer
IHC Recommend dilution:	30-150

