

## 兔抗 HSD17B7 多克隆抗体

中文名称： 兔抗 HSD17B7 多克隆抗体

英文名称： Anti-HSD17B7 rabbit polyclonal antibody

别名： PRAP; SDR37C1

相关类别： 一抗

抗原： HSD17B7

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	HSD17B7 encodes an enzyme that functions both as a 17-beta-hydroxysteroid dehydrogenase in the biosynthesis of sex steroids and as a 3-ketosteroid reductase in the biosynthesis of cholesterol.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	HSD17B7
<b>Immunogen:</b>	Fusion protein of human HSD17B7
<b>Full name:</b>	hydroxysteroid (17-beta) dehydrogenase 7
<b>Synonyms :</b>	PRAP; SDR37C1
<b>SwissProt:</b>	P56937
<b>ELISA Recommended dil</b>	1000-2000

<b>ution:</b>	
<b>WB Predicted band size:</b>	38 kDa
<b>WB Positive control:</b>	Human testis tissue
<b>WB Recommended diluti on:</b>	200-1000

