

兔抗 MT-ND4L 多克隆抗体

中文名称：兔抗 MT-ND4L 多克隆抗体

英文名称：Anti-MT-ND4L rabbit polyclonal antibody

别名：mitochondrially encoded NADH 4L dehydrogenase; MTND4L; ND4L

相关类别：一抗

储存：冷冻（-20℃）

抗原：MT-ND4L

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I) that is believed to belong to the minimal assembly required for catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone (By similarity).
Applications:	WB
Name of antibody:	MT-ND4L
Immunogen:	Synthetic peptide of human MT-ND4L
Full name:	mitochondrially encoded NADH 4L dehydrogenase
Synonyms:	MTND4L; ND4L

SwissProt:	P03901
WB Predicted band size:	11 kDa
WB Positive control:	MCF7 whole cell lysates
WB Recommended dilution:	500-1000

