

兔抗 MT-ND5 多克隆抗体

中文名称:	兔抗 MT-ND5 多克隆抗体
英文名称:	Anti-MT-ND5 rabbit polyclonal antibody
别名:	MTND5; ND5
储存:	冷冻 (-20℃)
抗原:	MT-ND5
宿主:	Rabbit
反应种属:	Human
相关类别:	一抗
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	Core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I) that is believed to belong to the minimal assembly required for catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.
Applications:	ELISA, IHC
Name of antibody:	MT-ND5
Immunogen:	Synthetic peptide of human MT-ND5
Full name:	mitochondrially encoded NADH dehydrogenase 5
Synonyms:	MTND5; ND5

SwissProt:	P03915
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and human ovarian cancer
IHC Recommend dilution:	25-100



