

## 兔抗 PTPN22 多克隆抗体

中文名称：兔抗 PTPN22 多克隆抗体

英文名称： Anti-PTPN22 rabbit polyclonal antibody

别名： LYP; PEP; LYP1; LYP2; PTPN8

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PTPN22

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene encodes of member of the non-receptor class 4 subfamily of the protein-tyrosine phosphatase family. The encoded protein is a lymphoid-specific intracellular phosphatase that associates with the molecular adapter protein CBL and may be involved in regulating CBL function in the T-cell receptor signaling pathway. Mutations in this gene may be associated with a range of autoimmune disorders including Type 1 Diabetes, rheumatoid arthritis, systemic lupus erythematosus and Graves' disease. Alternatively spliced transcript variants encoding distinct isoforms

	ms have been described.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PTPN22
<b>Immunogen:</b>	Fusion protein of human PTPN22
<b>Full name:</b>	protein tyrosine phosphatase, non-receptor type 22 (lymphoid)
<b>Synonyms :</b>	LYP; PEP; LYP1; LYP2; PTPN8
<b>SwissProt:</b>	Q9Y2R2
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and human ovarian cancer
<b>IHC Recommend dilution:</b>	30-150

