

兔抗 HSDL2 多克隆抗体

- 中文名称：兔抗 HSDL2 多克隆抗体
- 英文名称：Anti-HSDL2 rabbit polyclonal antibody
- 别名：hydroxysteroid dehydrogenase like 2; C9orf99; SDR13C1
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：HSDL2
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:	Has apparently no steroid dehydrogenase activity.
Applications:	ELISA, WB, IHC
Name of antibody:	HSDL2
Immunogen:	Fusion protein of human HSDL2
Full name:	hydroxysteroid dehydrogenase like 2
Synonyms:	C9orf99; SDR13C1
SwissProt:	Q6YN16
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	50-100
WB Predicted band size:	45 kDa

WB Positive control:	293T and Hela cell lysates
WB Recommended dilution:	1000-5000



