

兔抗 HSF1(Ab-142) 多克隆抗体

中文名称: 兔抗 HSF1(Ab-142) 多克隆抗体

英文名称: Anti-HSF1(Ab-142) rabbit polyclonal antibody

别名: HSTF1

储存: 冷冻 (-20℃) 避光

抗原: HSF1(Ab-142)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The product of this gene is a heat-shock transcription factor. Transcription of heat-shock genes is rapidly induced after temperature stress. Hsp90, by itself and/or associated with multichaperone complexes, is a major repressor of this gene.
Applications:	WB
Name of antibody:	HSF1(Ab-142)
Immunogen:	Synthesized non-phosphopeptide derived from human HSF1 around the phosphorylation site of threonine 142 (L-L-T(p)-D-V).
Full name:	heat shock transcription factor 1
Synonyms :	HSTF1
SwissProt:	Q00613

WB Predicted band size:	57 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000

