

## 兔抗 GDPD5 多克隆抗体

中文名称：兔抗 GDPD5 多克隆抗体

英文名称： Anti-GDPD5 rabbit polyclonal antibody

别名： GDE2; PP1665

相关类别： 一抗

储存： 冷冻（-20℃）

克隆类型： rabbit polyclonal

抗原： GDPD5

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

### 技术规格

<b>Background:</b>	Glycerophosphodiester phosphodiesterases, such as GDPD5 , are involved in glycerol metabolism. Promotes neurite for mation. Cooperates with PRDX1 to drive postmitotic motor neuron differentiation. The glycerophosphodiester phospho diesterase activity may be required for its role in neuronal differentiation. May contribute to the osmotic regulation o f cellular glycerophosphocholine.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	GDPD5
<b>Immunogen:</b>	Fusion protein of human GDPD5
<b>Full name:</b>	glycerophosphodiester phosphodiesterase domain containin

	g 5
<b>Synonyms :</b>	GDE2; PP1665
<b>SwissProt:</b>	Q8WTR4
<b>ELISA Recommended dilution :</b>	2000-5000
<b>WB Predicted band size:</b>	69 kDa
<b>WB Positive control:</b>	Mouse pancreas and human fetal brain tissue
<b>WB Recommended dilution:</b>	500-2000

