

兔抗 HSF2BP 多克隆抗体

中文名称: 兔抗 HSF2BP 多克隆抗体

英文名称: Anti-HSF2BP rabbit polyclonal antibody

别名: heat shock transcription factor 2 binding protein

储存: 冷冻 (-20℃)

抗原: HSF2BP

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	HSF2 binding protein (HSF2BP) associates with HSF2. The interaction occurs between the trimerization domain of HSF2 and the amino terminal hydrophilic region of HSF2BP that comprises two leucine zipper motifs. HSF2BP may therefore be involved in modulating HSF2 activation.
Applications:	ELISA, IHC
Name of antibody:	HSF2BP
Immunogen:	Synthetic peptide of human HSF2BP
Full name:	heat shock transcription factor 2 binding protein
SwissProt:	O75031
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human tonsil and Human ovarian cancer
IHC Recommend dilution:	30-150



