

兔抗 HSF2 多克隆抗体

中文名称：兔抗 HSF2 多克隆抗体

英文名称：Anti-HSF2 rabbit polyclonal antibody

别名：HSF 2, HSTF 2

抗原：HSF2

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the HSF family of transcription factors that bind specifically to the heat-shock promoter element and activate transcription. Heat shock transcription factors activate heat-shock response genes under conditions of heat or other stresses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB
Name of antibody:	HSF2
Immunogen:	Fusion protein of human HSF2
Full name:	heat shock transcription factor 2

Synonyms :	HSF 2, HSTF 2
SwissProt:	Q03933
ELISA Recommended dilution:	1000-5000
WB Predicted band size:	60 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	500-2000

