

兔抗 HSFX1 多克隆抗体

中文名称：兔抗 HSFX1 多克隆抗体

英文名称：Anti-HSFX1 rabbit polyclonal antibody

别名：LW-1

抗原：HSFX1

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Belongs to the HSF family.
Applications:	WB
Name of antibody:	HSFX1
Immunogen:	Synthesized peptide derived from N-terminal of human LW-1
Full name:	heat shock transcription factor family, X linked 1
Synonyms :	LW-1
SwissProt:	Q9UBD0
WB Predicted band size:	47 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution	500-3000

:

