

兔抗 GEMIN7 多克隆抗体

中文名称：兔抗 GEMIN7 多克隆抗体

英文名称： Anti-GEMIN7 rabbit polyclonal antibody

别名： gem nuclear organelle associated protein 7; SIP3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GEMIN7

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a component of the core SMN complex, which is required for pre-mRNA splicing in the nucleus. The encoded protein is found in the nucleoplasm, in nuclear "gems" (Gemini of Cajal bodies), and in the cytoplasm. Three transcript variants encoding the same protein have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	GEMIN7
Immunogen:	Fusion protein of human GEMIN7
Full name:	gem nuclear organelle associated protein 7
Synonyms:	SIP3

SwissProt:	Q9H840
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-200

