

兔抗 MTHFD2L 多克隆抗体

中文名称: 兔抗 MTHFD2L 多克隆抗体

英文名称: Anti-MTHFD2L rabbit polyclonal antibody

别名: methylenetetrahydrofolate dehydrogenase (NADP+ dependent) 2 like

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: MTHFD2L

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The specific function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	MTHFD2L
Immunogen:	Fusion protein of human MTHFD2L
Full name:	methylenetetrahydrofolate dehydrogenase (NADP+ dependent)) 2 like
SwissProt:	Q9H903
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human lung cancer
IHC Recommend dilution:	30-150



