

兔抗 PTPRM 多克隆抗体

- 中文名称：兔抗 PTPRM 多克隆抗体
- 英文名称：Anti-PTPRM rabbit polyclonal antibody
- 别名：RPTPM, RPTPU, PTPRL1, hR-PTPu, R-PTP-MU
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：PTPRM
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family . PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP possesses an extracellular region, a single transmembrane region, and two tandem catalytic domains, and thus represents a receptor-type PTP. The extracellular region contains a meprin-A5 antigen-PTP mu (

	<p>MAM) domain, an Ig-like domain and four fibronectin type III-like repeats. This PTP has been shown to mediate cell-cell aggregation through the interaction with another molecule of this PTP on an adjacent cell. This PTP can interact with scaffolding protein RACK1/GNB2L1, which may be necessary for the downstream signaling in response to cell-cell adhesion. Alternative splicing results in multiple transcripts encoding distinct isoforms.</p>
Applications:	ELISA, IHC
Name of antibody:	PTPRM
Immunogen:	Synthetic peptide of human PTPRM
Full name:	Protein tyrosine phosphatase, receptor type, M
Synonyms :	RPTPM, RPTPU, PTPRL1, hR-PTPu, R-PTP-MU
SwissProt:	P28827
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human breast cancer and Human lung cancer
IHC Recommend dilution:	50-200

