

DPYSL3 抗原（重组蛋白）

中文名称：DPYSL3 抗原（重组蛋白）

英文名称：DPYSL3 Antigen (Recombinant Protein)

别名：dihydropyrimidinase like 3; DRP3; ULIP; CRMP4; DRP-3; LCRMP; CRMP-4; ULIP-1

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 2-121 amino acids of human DPYSL3

技术规格

Full name:	dihydropyrimidinase like 3
Synonyms:	DRP3; ULIP; CRMP4; DRP-3; LCRMP; CRMP-4; ULIP-1
Swissprot:	Q14195
Gene Accession:	BC039006
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon guidance, neuronal growth cone collapse and cell migration (By similarity).