

## 兔抗 PTPRN 多克隆抗体

中文名称：兔抗 PTPRN 多克隆抗体

英文名称：Anti-PTPRN rabbit polyclonal antibody

别名：IA2; IA-2; ICA512; R-PTP-N; IA-2/PTP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PTPRN

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP possesses an extracellular region, a single transmembrane region, and a single catalytic domain, and thus represents a receptor-type PTP.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	PTPRN
<b>Immunogen:</b>	Fusion protein of human PTPRN
<b>Full name:</b>	protein tyrosine phosphatase, receptor type, N
<b>Synonyms :</b>	IA2; IA-2; ICA512; R-PTP-N; IA-2/PTP
<b>SwissProt:</b>	Q16849
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

