

兔抗 MTHFR 多克隆抗体

中文名称：兔抗 MTHFR 多克隆抗体

英文名称：Anti-MTHFR rabbit polyclonal antibody

储存：冷冻（-20℃）

抗原：MTHFR

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene catalyzes the conversion of 5,10-methylenetetrahydrofolate to 5-methyltetrahydrofolate, a co-substrate for homocysteine remethylation to methionine. Genetic variation in this gene influences susceptibility to occlusive vascular disease, neural tube defects, colon cancer and acute leukemia, and mutations in this gene are associated with methylenetetrahydrofolate reductase deficiency.
Applications:	ELISA, IHC
Name of antibody:	MTHFR
Immunogen:	Synthetic peptide of human MTHFR
Full name:	methylenetetrahydrofolate reductase (NAD(P)H)
SwissProt:	P42898

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human gastric cancer
IHC Recommend dilution:	25-100



