

兔抗 PTPRZ1 多克隆抗体

中文名称：兔抗 PTPRZ1 多克隆抗体

英文名称： Anti-PTPRZ1 rabbit polyclonal antibody

别名： protein tyrosine phosphatase, receptor-type, Z polypeptide 1

储存： 冷冻（-20℃） 避光

抗原： PTPRZ1

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the receptor protein tyrosine phosphatase family. Expression of this gene is restricted to the central nervous system (CNS), and it may be involved in the regulation of specific developmental processes in the CNS. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
Applications:	WB
Name of antibody:	PTPRZ1

Immunogen:	Synthesized peptide derived from N-terminal of human PTPRZ1.
Full name:	protein tyrosine phosphatase, receptor-type, Z polypeptide 1
Synonyms :	protein tyrosine phosphatase, receptor-type, Z polypeptide 1
SwissProt:	P23471
WB Predicted band size:	255 kDa
WB Positive control:	A549 cells lysate
WB Recommended dilution:	500-3000

