

## 兔抗 HSPA1L 多克隆抗体

中文名称: 兔抗 HSPA1L 多克隆抗体

英文名称: Anti-HSPA1L rabbit polyclonal antibody

别名: HSP70T; hum70t; HSP70-1L; HSP70-HOM

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: HSPA1L

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a 70kDa heat shock protein. In conjunction with other heat shock proteins, this protein stabilizes existing proteins against aggregation and mediates the folding of newly translated proteins in the cytosol and in organelles. The gene is located in the major histocompatibility complex class III region, in a cluster with two closely related genes which also encode isoforms of the 70kDa heat shock protein.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	HSPA1L

<b>Immunogen:</b>	Fusion protein of human HSPA1L
<b>Full name:</b>	heterogeneous nuclear ribonucleoprotein K
<b>Synonyms :</b>	HSP70T; hum70t; HSP70-1L; HSP70-HOM
<b>SwissProt:</b>	P34931
<b>WB Predicted band size:</b>	70 kDa
<b>WB Positive control:</b>	Hela, THP1 and ES-2 cells, testis tissue
<b>WB Recommended dilution:</b>	500-2000

