

兔抗 MTM1 多克隆抗体

中文名称：兔抗 MTM1 多克隆抗体

英文名称： Anti-MTM1 rabbit polyclonal antibody

别名： myotubularin 1; CNM; MTMX; XLMTM

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MTM1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a dual-specificity phosphatase that acts on both phosphotyrosine and phosphoserine. It is required for muscle cell differentiation and mutations in this gene have been identified as being responsible for X-linked myotubular myopathy.
Applications:	ELISA, WB
Name of antibody:	MTM1
Immunogen:	Fusion protein of human MTM1
Full name:	myotubularin 1
Synonyms:	CNM; MTMX; XLMTM
SwissProt:	Q13496

ELISA Recommended dilution:	5000-10000
WB Predicted band size:	70 kDa
WB Positive control:	Mouse thymus tissue lysate
WB Recommended dilution:	200-1000

