

兔抗 MTMR4 多克隆抗体

中文名称：兔抗 MTMR4 多克隆抗体

英文名称： Anti-MTMR4 rabbit polyclonal antibody

别名： ZFYVE11; FYVE-DSP2

相关类别： 一抗

储存： 冷冻 (-20℃)

抗原： MTMR4

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	MTMR4(myotubularin-related protein-4) cDNA and cDNAs encoding 3 other novel members of the myotubularin protein family. Dephosphorylates proteins phosphorylated on Ser, Thr, and Tyr residues and low molecular weight phosphatase substrate para-nitrophenylphosphate. Phosphorylates phosphatidylinositol 3,4,5-trisphosphate (PIP3).
Applications:	ELISA, IHC
Name of antibody:	MTMR4
Immunogen:	Fusion protein of human MTMR4
Full name:	myotubularin related protein 4
Synonyms :	ZFYVE11; FYVE-DSP2

SwissProt:	Q9NYA4
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	15-50

