

兔抗 PTTG1IP 多克隆抗体

中文名称：兔抗 PTTG1IP 多克隆抗体

英文名称： Anti-PTTG1IP rabbit polyclonal antibody

别名： PBF; C21orf1; C21orf3

相关类别： 一抗

抗原： PTTG1IP

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a single-pass type I integral membrane protein, which binds to pituitary tumor-transforming 1 protein (PTTG1), and facilitates translocation of PTTG1 into the nucleus. Coexpression of this protein and PTTG1 induces transcriptional activation of basic fibroblast growth factor. Alternatively spliced transcript variants have been found for this gene.
Applications:	WB, IHC

Name of antibody:	PTTG1IP
Immunogen:	Synthesized peptide derived from C-terminal of human PTTG.
Full name:	pituitary tumor-transforming 1 interacting protein
Synonyms :	PBF; C21orf1; C21orf3
SwissProt:	P53801
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	20 kDa
WB Positive control:	RAW264.7 cells lysate
WB Recommended dilution:	500-3000



