

兔抗 MTNR1B 多克隆抗体

中文名称：兔抗 MTNR1B 多克隆抗体

英文名称： Anti-MTNR1B rabbit polyclonal antibody

别名： MT2; FGQTL2; MEL-1B-R

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MTNR1B

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes one of two high affinity forms of a receptor for melatonin, the primary hormone secreted by the pineal gland. This gene product is an integral membrane protein that is a G-protein coupled, 7-transmembrane receptor. It is found primarily in the retina and brain although its detection requires RT-PCR. It is thought to participate in light-dependent functions in the retina and may be involved in the neurobiological effects of melatonin.
Applications:	ELISA, IHC
Name of antibody:	MTNR1B
Immunogen:	Synthetic peptide of human MTNR1B

Full name:	melatonin receptor 1B
Synonyms:	MT2; FGQTL2; MEL-1B-R
SwissProt:	P49286
ELISA Recommended dilution :	5000-10000
IHC positive control:	Human brain and human ovarian cancer
IHC Recommend dilution:	30-150



