

兔抗 PUM3 多克隆抗体

中文名称：兔抗 PUM3 多克隆抗体

英文名称：Anti-PUM3 rabbit polyclonal antibody

别名：pumilio RNA binding family member 3; PEN; HA-8; PUF6; XTP5; PUF-A; HLA-HA8; KIAA0020

抗原：PUM3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Inhibits the poly(ADP-ribosyl)ation activity of PARP1 and the degradation of PARP1 by CASP3 following genotoxic stress (PubMed:21266351). Binds to double-stranded RNA or DNA without sequence specificity (PubMed:25512524). Involved in development of the eye and of primordial germ cells (By similarity).
Applications:	ELISA, WB, IHC
Name of antibody:	PUM3

Immunogen:	Synthetic peptide of human PUM3
Full name:	pumilio RNA binding family member 3
Synonyms:	PEN; HA-8; PUF6; XTP5; PUF-A; HLA-HA8; KIAA0020
SwissProt:	Q15397
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	74 kDa
WB Positive control:	Raji cell lysate
WB Recommended dilution:	1000-5000



