

## 兔抗 GGCX 多克隆抗体

中文名称：兔抗 GGCX 多克隆抗体

英文名称： Anti-GGCX rabbit polyclonal antibody

别名： VKCFD1

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： GGCX

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes an enzyme which catalyzes the posttranslational modification of vitamin K-dependent protein. Many of these vitamin K-dependent proteins are involved in coagulation so the function of the encoded enzyme is essential for hemostasis. Mutations in this gene are associated with vitamin K-dependent coagulation defect and PXE-like disorder with multiple coagulation factor deficiency. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GGCX

<b>Immunogen:</b>	Fusion protein of human GGCX
<b>Full name:</b>	gamma-glutamyl carboxylase
<b>Synonyms :</b>	VKCFD1
<b>SwissProt:</b>	P38435
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human breast cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	100-300

