

兔抗 MTOR 多克隆抗体

中文名称：兔抗 MTOR 多克隆抗体

英文名称： Anti-MTOR rabbit polyclonal antibody

别名： FRAP, FRAP1, FRAP2, RAFT1, RAPT1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MTOR

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to a family of phosphatidylinositol kinase-related kinases. These kinases mediate cellular responses to stresses such as DNA damage and nutrient deprivation. This protein acts as the target for the cell-cycle arrest and immunosuppressive effects of the FKBP12-rapamycin complex. The ANGPTL7 gene is located in an intron of this gene.
Applications:	ELISA, IHC
Name of antibody:	MTOR
Immunogen:	Synthetic peptide of human MTOR
Full name:	mechanistic target of rapamycin (serine/threonine kinase)

Synonyms :	FRAP, FRAP1, FRAP2, RAFT1, RAPT1
SwissProt:	P42345
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

