

兔抗 PWWP2B 多克隆抗体

中文名称: 兔抗 PWWP2B 多克隆抗体

英文名称: Anti-PWWP2B rabbit polyclonal antibody

别名: PWWP2; pp8607; bA432J24.1

储存: 冷冻 (-20℃) 避光

抗原: PWWP2B

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	a protein of unknown function.
Applications:	WB
Name of antibody:	PWWP2B
Immunogen:	Synthesized peptide derived from internal of human PWWP2B.
Full name:	PWWP domain containing 2B
Synonyms :	PWWP2; pp8607; bA432J24.1
SwissProt:	Q6NUJ5
WB Predicted band size:	64 kDa

WB Positive control:	LOVO cells lysate
WB Recommended dilution:	500-3000

