

兔抗 MTTP 多克隆抗体

中文名称：兔抗 MTTP 多克隆抗体

英文名称：Anti-MTTP rabbit polyclonal antibody

别名：ABL; MTP

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：MTTP

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	MTP encodes the large subunit of the heterodimeric microsomal triglyceride transfer protein. Protein disulfide isomerase (PDI) completes the heterodimeric microsomal triglyceride transfer protein, which has been shown to play a central role in lipoprotein assembly. Mutations in MTP can cause abetalipoproteinemia.
Applications:	WB
Name of antibody:	MTTP
Immunogen:	Fusion protein of human MTTP
Full name:	microsomal triglyceride transfer protein
Synonyms :	ABL; MTP

SwissProt:	P55157
WB Predicted band size:	99 kDa
WB Positive control:	293 cells
WB Recommended dilution:	500-2000

