

兔抗 MTUS1 多克隆抗体

中文名称：兔抗 MTUS1 多克隆抗体

英文名称： Anti-MTUS1 rabbit polyclonal antibody

别名： ATBP; ATIP; ICIS; MP44; MTSG1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MTUS1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein which contains a C-terminal domain able to interact with the angiotension II (AT2) receptor and a large coiled-coil region allowing dimerization. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. One of the transcript variants has been shown to encode a mitochondrial protein that acts as a tumor suppressor and participates in AT2 signaling pathways.
Applications:	ELISA, IHC
Name of antibody:	MTUS1
Immunogen:	Fusion protein of human MTUS1

Full name:	microtubule associated tumor suppressor 1
Synonyms :	ATBP; ATIP; ICIS; MP44; MTSG1
SwissProt:	Q9ULD2
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human breast cancer
IHC Recommend dilution:	50-200

