

## 兔抗 GHDC 多克隆抗体

中文名称：兔抗 GHDC 多克隆抗体

英文名称： Anti-GHDC rabbit polyclonal antibody

别名： GH3 domain containing; LGP1; D11LGP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GHDC

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	GHDC, also known as LGP1, is a GH3 domain-containing protein that is highly expressed in normal mammary tissues and mammary tumors. GHDC is primarily located in the nuclear envelope. In mouse, its gene is located on chromosome 11 adjacent to the STAT3/5 locus. Little is known of the function of this protein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GHDC
<b>Immunogen:</b>	Synthetic peptide of human GHDC
<b>Full name:</b>	GH3 domain containing
<b>Synonyms:</b>	LGP1; D11LGP1

<b>SwissProt:</b>	Q8N2G8
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-300



