

兔抗 PXMP2 多克隆抗体

中文名称: 兔抗 PXMP2 多克隆抗体

英文名称: Anti-PXMP2 rabbit polyclonal antibody

别名: PMP22

储存: 冷冻 (-20℃) 避光

抗原: PXMP2

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane.
Applications:	IHC
Name of antibody:	PXMP2
Immunogen:	Synthesized peptide derived from internal of human PXMP2.
Full name:	peroxisomal membrane protein 2, 22kDa
Synonyms :	PMP22

SwissProt:	Q9NR77
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

