

## 兔抗 PXN (Phospho-Tyr118)多克隆抗体

- 中文名称： 兔抗 PXN (Phospho-Tyr118)多克隆抗体
- 英文名称： Anti-PXN (Phospho-Tyr118) rabbit polyclonal antibody
- 相关类别： 一抗
- 储 存： 冷冻（-20℃） 避光
- 宿 主： Rabbit
- 抗 原： PXN (Phospho-Tyr118)
- 反应种属： Human, Mouse, Rat
- 标 记 物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	PXN (Phospho-Tyr118)
<b>Immunogen:</b>	Synthetic peptide of human PXN (Phospho-Tyr118)
<b>Full name:</b>	paxillin (Phospho-Tyr118)
<b>SwissProt:</b>	P49023

<b>IHC positive control:</b>	Human breast carcinoma
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	65 kDa
<b>WB Positive control:</b>	NIH/3T3 cells untreated or treated with serum
<b>WB Recommended dilution:</b>	500-1000



