

兔抗 PXN(Ab-272) 多克隆抗体

中文名称：兔抗 PXN(Ab-272) 多克隆抗体

英文名称： Anti-PXN(Ab-272) rabbit polyclonal antibody

抗 原： PXN(Ab-272)

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extra cellular matrix (focal adhesion).
Applications:	WB
Name of antibody:	PXN(Ab-272)
Immunogen:	Synthesized non-phosphopeptide derived from human Paxillin around the phosphorylation site of serine 272 (M-A-S(p)-L-S).
Full name:	paxillin
Synonyms :	0
SwissProt:	P49023

WB Predicted band size:	65 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000

