

## 兔抗 HSPBAP1 多克隆抗体

中文名称：兔抗 HSPBAP1 多克隆抗体

英文名称：Anti-HSPBAP1 rabbit polyclonal antibody

别名：PASS1

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：HSPBAP1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a protein that binds to one of the small heat shock proteins, specifically hsp27. Hsp27 is involved with cell growth and differentiation. This encoded protein was found to be abnormally expressed in patients with intractable epilepsy, although how brain function is affected remains unknown.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	HSPBAP1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human HBAP1.
<b>Full name:</b>	HSPB (heat shock 27kDa) associated protein 1
<b>Synonyms :</b>	PASS1

<b>SwissProt:</b>	Q96EW2
<b>WB Predicted band size:</b>	55 kDa
<b>WB Positive control:</b>	Hela cells and 293 cells lysates
<b>WB Recommended dilution:</b>	500-3000

