

兔抗 MUC15 多克隆抗体

中文名称：兔抗 MUC15 多克隆抗体

英文名称： Anti-MUC15 rabbit polyclonal antibody

别名： PAS3; MUC-15; PASIII

抗原： MUC15

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Mucins are a family of high molecular weight, heavily glycosylated proteins (glycoconjugates) produced by epithelial tissues in most metazoans. Mucin 15 is involved in cell adhesion to the extracellular matrix. It is primarily expressed in prostate, testis, lung, spleen, thymus, ovary, small intestine, colon, bone marrow, lymph node, lung and peripheral blood leukocytes. Mucin 15 is a highly glycosylated protein.
Applications:	ELISA, IHC
Name of antibody:	MUC15
Immunogen:	Fusion protein of human MUC15
Full name:	mucin 15, cell surface associated

Synonyms :	PAS3; MUC-15; PASIII
SwissProt:	Q8N387
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human lung cancer
IHC Recommend dilution:	25-100

