

兔抗 MUC16 多克隆抗体

中文名称：兔抗 MUC16 多克隆抗体

英文名称：Anti-MUC16 rabbit polyclonal antibody

别名：CA125

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MUC16

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Synthetic peptide of human MUC16
Applications:	ELISA, IHC
Name of antibody:	MUC16
Immunogen:	Synthetic peptide of human MUC16
Full name:	mucin 16, cell surface associated
Synonyms:	CA125
SwissProt:	Q8WXI7
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer;
IHC Recommend dilution:	15-50

