

兔抗 MUC2 多克隆抗体

中文名称：兔抗 MUC2 多克隆抗体

英文名称：Anti-MUC2 rabbit polyclonal antibody

别名：MLP, SMUC, MUC-2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MUC2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the mucin protein family. Mucins are high molecular weight glycoproteins produced by many epithelial tissues. The protein encoded by this gene is secreted and forms an insoluble mucous barrier that protects the gut lumen. The protein polymerizes into a gel of which 80% is composed of oligosaccharide side chains by weight. The protein features a central domain containing tandem repeats rich in threonine and proline that varies between 50 and 115 copies in different individuals. Alternatively spliced transcript variants of this gene have been described, but their full-length nature is

	not known.
Applications:	ELISA, IHC
Name of antibody:	MUC2
Immunogen:	Synthetic peptide of human MUC2
Full name:	mucin 2, oligomeric mucus/gel-forming
Synonyms :	MLP, SMUC, MUC-2
SwissProt:	Q02817
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer and Human lung cancer
IHC Recommend dilution:	50-200

