

## 兔抗 MUC3A 多克隆抗体

中文名称：兔抗 MUC3A 多克隆抗体

英文名称： Anti-MUC3A rabbit polyclonal antibody

别名： MUC3, MUC-3A

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MUC3A

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                          |   |
|--------------------------|---|
| <b>Background:</b>       | The mucin genes encode epithelial glycoproteins, some of which are secreted and some membrane bound. Each of the genes contains at least one large domain of tandemly repeated sequence that encodes the peptide sequence rich in serine and/or threonine residues, which carries most of the O-linked glycosylation. |
| <b>Applications:</b>     | ELISA, IHC  |
| <b>Name of antibody:</b> | MUC3A   |
| <b>Immunogen:</b>        | Synthetic peptide of human MUC3A  |
| <b>Full name:</b>        | mucin 3A, cell surface associated   |
| <b>Synonyms :</b>        | MUC3, MUC-3A  |

|                                    |  |
|------------------------------------|--|
| <b>SwissProt:</b>                  | Q02505                                   |
| <b>ELISA Recommended dilution:</b> | 1000-5000                                |
| <b>IHC positive control:</b>       | Human colon cancer and Human lung cancer |
| <b>IHC Recommend dilution:</b>     | 25-100                                   |

