

兔抗 MUC5AC 多克隆抗体

中文名称：兔抗 MUC5AC 多克隆抗体

英文名称： Anti-MUC5AC rabbit polyclonal antibody

别名： TBM, leB, MUC5

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MUC5AC

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Mucin-5AC is a protein that in humans is encoded by the MUC5AC gene. This gene has been linked to mucous hypersecretion in the pulmonary tracts and associated to chronic obstructive pulmonary disease (COPD). Gel-forming glycoprotein of gastric and respiratory tract epithelia that protects the mucosa from infection and chemical damage by binding to inhaled microorganisms and particles that are subsequently removed by the mucociliary system. Highly expressed in surface mucosal cells of respiratory tract and stomach epithelia. Overexpressed in a number of carcinomas. Also expressed in Barrett's esophagus epithelium

	and in the proximal duodenum.
Applications:	ELISA, IHC
Name of antibody:	MUC5AC
Immunogen:	Synthetic peptide of human MUC5AC
Full name:	mucin 5AC, oligomeric mucus/gel-forming
Synonyms :	TBM, IeB, MUC5
SwissProt:	P98088
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	50-200

