

## 兔抗 MUSK 多克隆抗体

中文名称：兔抗 MUSK 多克隆抗体

英文名称： Anti-MUSK rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

抗 原： MUSK

宿 主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a muscle-specific tyrosine kinase receptor. The encoded protein may play a role in clustering of the acetylcholine receptor in the postsynaptic neuromuscular junction. Mutations in this gene have been associated with congenital myasthenic syndrome. Alternatively spliced transcript variants have been described.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MUSK
<b>Immunogen:</b>	Fusion protein of human MUSK
<b>Full name:</b>	muscle, skeletal, receptor tyrosine kinase
<b>SwissProt:</b>	O15146
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human liver cancer

IHC Recommend dilution: 50-200

