

兔抗 MVK 多克隆抗体

中文名称：兔抗 MVK 多克隆抗体

英文名称： Anti-MVK rabbit polyclonal antibody

别名： MK; LRBP; MVLK; POROK3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MVK

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes the peroxisomal enzyme mevalonate kinase. Mevalonate is a key intermediate, and mevalonate kinase a key early enzyme, in isoprenoid and sterol synthesis. Mevalonate kinase deficiency caused by mutation of this gene results in mevalonic aciduria, a disease characterized by psychomotor retardation, failure to thrive, hepatosplenomegaly, anemia and recurrent febrile crises. Defects in this gene also cause hyperimmunoglobulinemia D and periodic fever syndrome, a disorder characterized by recurrent episodes of fever associated with lymphadenopathy, arthralgia, gastrointestinal distress and skin rash.

Applications:	ELISA, WB, IHC
Name of antibody:	MVK
Immunogen:	Fusion protein of human MVK
Full name:	mevalonate kinase
Synonyms :	MK; LRBP; MVLK; POROK3
SwissProt:	Q03426
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human esophagus cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	42 kDa
WB Positive control:	Raji cells
WB Recommended dilution:	200-1000



