

兔抗 MYBL1 多克隆抗体

中文名称：兔抗 MYBL1 多克隆抗体

英文名称： Anti-MYBL1 rabbit polyclonal antibody

别名： AMYB; A-MYB

相关类别： 一抗

抗原： MYBL1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Strong transcriptional activator; DNA-binding protein that specifically recognize the sequence 5'-YAAC[GT]G-3'. Could have a role in the proliferation and/or differentiation of neurogenic, spermatogenic and B-lymphoid cells.
Applications:	WB
Name of antibody:	MYBL1
Immunogen:	Synthesized peptide derived from internal of human MYB-A.
Full name:	v-myb avian myeloblastosis viral oncogene homolog-like 1
Synonyms :	AMYB; A-MYB
SwissProt:	P10243

WB Predicted band size:	86 kDa
WB Positive control:	LOVO cells lysate
WB Recommended dilution:	500-3000

