

兔抗 MYB 多克隆抗体

中文名称：兔抗 MYB 多克隆抗体

英文名称：Anti-MYB rabbit polyclonal antibody

别名：efg; Cmyb; c-myb; c-myb_CDS

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MYB

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a transcription factor that is a member of the MYB family of transcription factor genes. The protein contains three domains, an N-terminal DNA-binding domain, a central transcriptional activation domain and a C-terminal domain involved in transcriptional repression. This protein plays an essential role in the regulation of hematopoiesis and may play a role in tumorigenesis. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB
Name of antibody:	MYB
Immunogen:	Synthetic peptide of human MYB

Full name:	v-myb avian myeloblastosis viral oncogene homolog
Synonyms :	efg; Cmyb; c-myb; c-myb_CDS
SwissProt:	P10242
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	72 kDa
WB Positive control:	Hela and Raji cells
WB Recommended dilution:	200-1000

