

兔抗 MYEF2 多克隆抗体

中文名称：兔抗 MYEF2 多克隆抗体

英文名称：Anti-MYEF2 rabbit polyclonal antibody

别名：myelin expression factor 2; MEF-2; MST156; myEF-2; MSTP156; HsT18564

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MYEF2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | Transcriptional repressor of the myelin basic protein gene (MBP). Binds to the proximal MB1 element 5'-TTGTCC-3' of the MBP promoter. Its binding to MB1 and function are inhibited by PURA (By similarity). |
| Applications: | ELISA, WB, IHC |
| Name of antibody: | MYEF2 |
| Immunogen: | Fusion protein of human MYEF2 |
| Full name: | myelin expression factor 2 |
| Synonyms: | MEF-2; MST156; myEF-2; MSTP156; HsT18564 |
| SwissProt: | Q9P2K5 |
| ELISA Recommended dilution: | 5000-10000 |

| | |
|---------------------------------|---|
| IHC positive control: | Human ovarian cancer and Human thyroid cancer |
| IHC Recommend dilution: | 50-200 |
| WB Predicted band size: | 64 kDa |
| WB Positive control: | 231 cell lysate |
| WB Recommended dilution: | 500-2000 |



