

兔抗 MYF5 多克隆抗体

中文名称：兔抗 MYF5 多克隆抗体

英文名称： Anti-MYF5 rabbit polyclonal antibody

别名： bHLHc2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MYF5

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Myogenic factor 5 is a protein that in humans is encoded by the MYF5 gene. It is a protein with a key role in regulating muscle differentiation or myogenesis. Without Myf5 and MyoD, myogenic cells will fail to progress normally during the determination stage of myogenesis. Myf5 belongs to a family of proteins known as myogenic regulatory factors (MRFs). These bHLH (basic helix loop helix) transcription factors act sequentially in myogenic differentiation. MRF family members include Myf5, MyoD (Myf3), myogenin, and MRF4 (Myf6).
Applications:	ELISA, IHC

Name of antibody:	MYF5
Immunogen:	Synthetic peptide of human MYF5
Full name:	myogenic factor 5
Synonyms :	bHLHc2
SwissProt:	P13349
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colon cancer and Human thyroid cancer
IHC Recommend dilution:	25-100

