

兔抗 MYF6 多克隆抗体

中文名称：兔抗 MYF6 多克隆抗体

英文名称：Anti-MYF6 rabbit polyclonal antibody

别名：CNM3; MRF4; myf-6; bHLHc4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MYF6

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a probable basic helix-loop-helix (bHLH) DNA binding protein involved in muscle differentiation. The encoded protein likely acts as a heterodimer with another bHLH protein. Defects in this gene are a cause of autosomal dominant centronuclear myopathy (ADCNM).
Applications:	ELISA, IHC
Name of antibody:	MYF6
Immunogen:	Synthetic peptide of human MYF6
Full name:	myogenic factor 6
Synonyms:	CNM3; MRF4; myf-6; bHLHc4

SwissProt:	P23409
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

